

A-Scan Plus Connect Recommended Computer Specs

1. **Processor**
Intel Core 2 Duo 2.20 GHz, Core i3 2.2 GHz or higher

2. **Memory**
3.0 GB, DDR2-667 SDRAM or higher

3. **Screen Resolution of Monitor**
High resolution monitor (minimum 1366 X 768 pixels;
recommended 1600 X 900 pixels or higher).

4. **Hard Disk**
Minimum 80 GB, 7200 RPM, recommended 160 GB or higher for
high patient flow site (more than 30 B-Scans per month).

5. **USB Ports**
Minimum of 2 USB 2.0

6. **Medical Grade Power Supply**

7. **Software and Operation System (OS)**
Windows 7[®], 8[®], Microsoft Windows[®] Vista or Microsoft
Windows[®] XP Professional Edition, SP2