

Keeler
3222 Phoenixville Pike, Bldg. 50
Malvern, PA 19355
800.523.5620 | 610.353.4350


Haag Streit Fitted Slit Lamp Breath Shield

Ref: AX17415



With the current global health crisis, there has been a huge demand from healthcare professionals for products to increase their protection from airborne pathogens. As such Keeler is introducing a “Haag Streit Fitted Breath Shield”, that will introduce extra protection between the clinician and patient when using a slit lamp.

The breath shield will fit any slit lamp, being 15.25” wide by 19” high gives great protection while still allowing easy access with a Volk lens. The shield is made of 1.5mm clear plastic - making it light but also sturdy.




Haag Streit Fitted Breath Shield

Part Number: AX17415

Standard Breath Shield

Part Number: EP39-80242



Haag Streit Fitted Breath Shield Fitting

The shield is designed to fit any slit lamp, including Haag Streit. The ocular cut out sits at the base of the oculars and allows for full pupillary distance adjustments.

