


# Keeler is helping you and your patients

Keeler is working hard to provide the products needed to keep both you and your patients protected and safe. Availability is continuously changing, so please call your sales representative or customer service at 800.523.5620 or email [customerservice@keelerusa.com](mailto:customerservice@keelerusa.com) for more information.

## Face Masks



### KN95 Masks\*

Item #: AX17399 10 per package

- Virus protection, anti-dust, smoke, gas, allergies, germs and PPE for adults
- FDA and CE certified manufacturer
- Non-woven polypropylene

**\*Image may not represent the actual mask you receive.**



### 3-Layer Medical Mask

Item #: AX17400 50 per package

- Anti-dust, smoke, gas, allergies, germs and PPE for adults
- FDA cleared
- Disposable material
- Effectively isolating bacteria, pollen, dust



### 3-Ply Youth Mask

10 per package

- Kids 3-layer mask with ear loops
- Anti-fog mask with 99.8% bacterial filtration efficiency.

AX17496B Blue mask

AX17496P Pink mask



### FC Clean Reusable Protective Mask

Item #: AX17438

- The outer layer is made of cleanroom compatible fabric with droplet resistance
- The inner layer is made of antibacterial fabric
- Long antibacterial effect even after washing it more than 50 times
- 5 year shelf life
- Non-medical
- Measures 9cm x 18cm

## Hand Sanitizers



### Soho Fresh Hand Sanitizer

Item #: AX17445

- 33.81oz (1 L) squeeze top bottle
- Contains 70% alcohol
- Kills 99.99% of germs



### GSL Organics Hand Sanitizer

Item #: AX17404

- 6oz. gel with pump
- Concentrated with 80% grain alcohol
- Non-drying and moisturizing blend with powerful essential oils, glycerin, aloe vera, and grain alcohol to kill germs and leave your hands feeling moisturized



### Instant Gel Hand Sanitizer

Item #: AX17450

- 16oz. pump bottle
- Contains 75% alcohol
- Kills 99.99% of germs


---

**Availability is limited and on a first come, first serve basis.**

**Due to high demand, these items are non-returnable.**

---

## Breath Shields



### X-Large Slit Lamp Breath Shield

Item #: EP39-80753

- Provides both clinician and patient protection while still allowing easy use of the slit lamp
- Goes between the binoculars and magnification block
- Fits Keeler slit lamps and any slit lamp with a 70mm diameter dovetail joint
- Made of 1.5mm thick clear plastic
- Measures 35cm x 30cm (13.7" x 11.8"), the perfect size for reaching around



### Large Fitted Slit Lamp Breath Shield

Item #: AX17415


- Provides both clinician and patient protection while still allowing easy use of the slit lamp
- The shield goes around the oculars
- Fits all slit lamps
- Made of 1.5mm thick clear plastic
- Measures 7.5"H x 6"W
- Ocular cut out allows for full pupillary distance adjustments



### Keeler BIO Face Shield

Item #: AX17448

- Scratch and defect free
- Fits Keeler All pupil II and Vantage Plus
- Shield is held in place by the snaps on the headband foam liner
- Made of 0.5mm-0.8mm thick clear polyethylene terephthalate
- Measures 275mm W x 190mm H (10.8" W x 7.5" H")



### The Phoropter Shield™

Item #: AX17447

- Designed to easily slip onto most A style phoropters
- Does not interfere with the use of the phoropter, or touching the nose of the patient, even when the PD is adjusted
- Made of 30mil polycarbonate composition allowing for easy installation and removal
- Measures approximately 7.5" x 7.5"

## Glasses



### Virtua™ Scratch-Resistant Safety Glasses

Item #: AX17437

- Clear and scratch-resistant lens
- Frameless, wraparound frames
- Universal size



### Ergonomic PVC Anti Fog Safety Goggles

Item #: AX17449

- Clear safety polycarbonate lens
- Anti Fog and saliva droplet blockers
- UV400
- Adjustable comfort strap

## Ear Savers



### Ear Saver Mask Adjuster

Item #: AX17478 10 per package

- The Ear Saver allows the user to adjust the tightness and fit of the ear loops on the mask to their desired fit and obtain a proper seal
- Fits any mask with elastic earloops

## Face Shields



### Face Shield

Item #: AX17416

- 12.387"W x 10.5"H clear shield covers the whole face
- Has 10.6"W x 1"H headband

**For more information, please contact your sales representative at 800.523.5620.**

## Disinfectants



### Betco Fight Bac RTU

- This quaternary-based, anti-bacterial cleaner
- Cuts through grease, heavy soil, and filth with quick drying times.

AX17433 32oz. Bottle  
AX17433J Trigger spray with straw



## Eye Patches



### MASK-IT Eye Patch

Item #: AX17024 500 per package

- Disposable and adjustable eye patch
- Pre-slit patch creates various cone depths for eye lashes

## Thermometers



### Non-Contact Infrared Digital Thermometer

Item #: AX17396

- In as little as 1-3 seconds, you can get an accurate temperature reading through the activation
- 2 temperature units for conversion - Fahrenheit or Celsius
- Long service life. Can be used more than 100,000 times with two AA batteries (not included)

## Gowns



### Level 2 SMS Gowns

- Medical level 2 SMS gowns
- Single-use
- Made from 3- layer SMS fabric.

AX17489 Small  
AX17490 Medium  
AX17491 Large  
AX17492 X-Large

For more information, please contact your sales representative at 800.523.5620.